


INSTRUCTIONS FW01—FINGER EXTENSION

INDICATIONS—Increased PIP Joint Extension

Contraindicated for fixed contractures—Patients with fragile skin—Severe edema

APPLICATION—With curved lateral wires proximal, slide finger into splint. Position distal finger band at or proximal to DIP joint. Dorsal pad proximal to PIP joint. Proximal band at metacarpal head.


TO FACILITATE DONNING—Gently squeezing ends of splint towards each other slightly, will open space between pads and allows finger to slide into proper position.(A)



ADJUST EXTENSION FORCE W/SPLINT REMOVED (B)

INCREASE—Bend ends away from each other

DECREASE—Bend ends towards each other



CLEANING—Hand wash with damp towel, cold water and mild soap. Thoroughly rinse and air dry completely before reapplying.

TO BE FIT BY A QUALIFIED HEALTHCARE PROFESSIONAL

Qualified Healthcare Professional to provide use and wearing instructions/guidelines.

PRECAUTIONS—Immediately contact your healthcare professional if you are experiencing pain, swelling, numbness, skin discoloration, or have any questions or concerns.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use