


FW02 Instructions

PIP SPRING EXTENSION—CAPENER DESIGN

Coiled springs provide dynamic extension while allowing active flexion



APPLICATION—Loosen strap—Slide finger into splint—Position proximal pad under metacarpal head—Align coiled springs with P.I.P. joint—Position distal band and strap proximal to D.I.P. joint—Extension force can be adjusted by tighten or loosening strap—Secure strap around finger and remove excess.

TO BE FIT BY A QUALIFIED HEALTHCARE PROFESSIONAL

Qualified Healthcare Professional to provide use and wearing instructions/guidelines/precautions RX.

Immediately contact your healthcare professional if you are experiencing pain, swelling, numbness, skin discoloration, or have any questions or concerns.

CLEANING—Hand wash with damp towel, cold water and mild soap. Thoroughly rinse and air dry completely before reapplying.

Contraindicated for fixed or boney contractures—Severe edema—Fragile skin—Etc.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use