


INSTRUCTIONS #10 Composite Spring to extend distal joints

Open straps. Position splint on top of hand with M.P. bar between M.P. and P.I.P. joints and hand band behind M.P. joints. Secure wrist strap and position pad. Secure hand strap. Secure finger strap. Put finger loops on fingers. Splint is designed to dynamically extend distal joints while allowing active flexion. Evaluate fit, function and instruct patient on use.

Modifications: Outrigger 1 & 4 have optional position for finger loops. Rubber bands can be added or subtracted. Clock spring outriggers and bands can be carefully bent to improve fit and function. Remove finger loops before removing splint. Clean with mild soap and damp towel. Air dry completely before reapplying.

Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritations or have any questions or concerns.



BUNNELL #10 Composite Spring to Extend Distal Joints

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use