


## INSTRUCTIONS #1-A Cock-up Splint with Outrigger

Apply outrigger to splint. With the outrigger angled down firmly insert wire ends into bracket on bottom side of splint. Loosen wrist strap. Open forearm strap. Straps can be removed from splint and repositioned to facilitate closure with right or left hand. Slide hand into splint and under wrist strap until comfortably positioned against hand piece. Hand piece can be angled to fit right or left hand. Side tabs can be reshaped. Position strap pad over wrist. Cross strap over top of wrist and secure. Secure forearm strap. Apply finger loops and attach rubber bands to hook on bottom of splint. Evaluate fit, function and instruct patient on use. Rubber bands can be added or subtracted and clock spring can be bent to modify fit and function of splint.

Clean with mild soap and damp towel. Air dry completely before reapplying.

Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritations or have any questions or concerns.



Bunnell #1A Cock-up with outrigger used to dorsiflex wrist & flex M.P. joints

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use