


INSTRUCTIONS #2B Knuckle Bender w/extra segment to flex distal joints

Open strap on second segment. Applying splint with rubber band attached only to hook on first and second segment (no force to flex joints) will facilitate proper positioning of splint on hand. With first and second segment facing forward slide hand into splint. Position palmar bar in arch of palm and first and second segment in proper position to flex joints.

Secure strap on second segment. Attach rubber bands to hook on hand section. Splint is designed to dynamically flex M.P. and distal joints while allowing active extension.

Evaluate fit, function and instruct patient on use. Rubber bands can be added or subtracted and the wires and bands can be carefully bent to modify the fit and function of the splint. Clean with mild soap and damp towel. Air dry completely before reapplying.

Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritations or have any questions or concerns.



Bunnell #2B Knuckle Bender Splint to flex the metacarpophalangeal joints with extra segment to flex the distal joints.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use