


INSTRUCTIONS #5 Finger Extension w/ Side Tabs & Clock Spring

Open strap. Slide finger into splint. Wrap strap around wrist and comfortably secure. Evaluate fit, function and instruct patient on use. Splint can be carefully bent to modify fit and adjust force.

Clean with mild soap and damp towel. Air dry completely before reapplying.

Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritations or have any questions or concerns.



BUNNELL #5 Finger Extension with Clock Spring to extend a finger or thumb.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use