


## INSTRUCTIONS #9A Reverse knuckle bender with outrigger

Attach finger loops to openings at top of hand band. Sizes medium, large and extra large have four attachment openings (one for each finger). Extra small and small have two openings requiring two finger loops per opening. Open wrist strap and slide hand into splint. Position M.P. extension bar under fingers between M.P. and P.I.P. joints. Position round pad in palm. Bring strap over top of wrist and secure. Finish attaching rubber band to hooks. Pass finger loops over top of outrigger and position on fingers. Splint is designed to extend M.P. joints and distal joints while allowing active flexion.

Evaluate fit, function and instruct patient on use.

Modifications: Remove from patient. Rubber bands can be added or subtracted. Wire and band can be carefully bent to improve fit and function of splint. Clean with mild soap and damp towel. Air dry completely before reapplying. Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritations or have any questions or concerns.



BUNNELL #9A Reverse Knuckle Bender Splint to extend the metacarpophalangeal joints with outrigger to extend the distal interphalangeal joint.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use